

# Considerations on *ex vivo* Conversion of Prodrugs during Bioanalysis

**Xiaohan Cai, Ph.D.**

Senior Staff Fellow

Division of Generic Drugs Study Integrity

Office of Study Integrity and Surveillance

CDER | US FDA

# Disclaimer




The views and opinions presented here represent those of the speaker and should not be considered to represent advice or guidance on behalf of the U.S. Food and Drug Administration


# Prodrugs and Common Types

- Prodrugs are designed to transform to active drug moieties *in vivo*
  - Address limitations of the active form in ADMET
  - Through a chemical/enzymatic process or combination
- Common prodrug classes:
  - Ester prodrugs – e.g. Fesoterodine
  - Phosphates prodrugs – e.g. Tedizolid phosphate
  - Amides prodrugs – e.g. Tenofovir alafenamide
  - Carbamate prodrugs – e.g. Loratadine

# Prodrugs in Biological Matrices



# What happen to Study Samples *ex vivo*?



## Considerations during Method Development

- Understand the chemical and/or enzymatic pathways leading to the conversion
  - Different stability profiles in varied species and matrices
- Rigorous stability testing
  - In applicable biological matrices
  - Impact from sample extractions
    - Processing temperature
    - Type of extractions: PPT, LLE, SPE, etc
    - Extraction solutions
  - Different storage conditions: benchtop, frozen, F/T cycles, post-extraction
- Possible in-source fragmentation of the prodrug when using an LC-MS/MS method

## Considerations during Method Development

- If possible, monitoring the changes for both the prodrug and active drug under stressed conditions
  - Provides information on the rate and level of conversion
- Consider the range of concentration ratio between the prodrug and the active drug in the in vivo study
  - Critical information to decide the possible data impact for the in vivo studies
  - For example, 10% *ex vivo* conversion has different impact on concentration change for the active drug when:
    - Sample with prodrug/active drug molar ratio at 1:1
    - Sample with prodrug/active drug molar ratio at 10:1

## Things to Consider during Method Development

- When needed, test different approaches to stabilize prodrugs in biological samples ex vivo. Common approaches include:
  - Temperature control
  - pH control
  - Enzyme inhibitors
  - Vacutainer tubes containing inhibitors
  - Other effective approach




## Consideration during Method Validation and Sample Collection

- Method Validation: Parameters validated according to the latest BMV guidance. Valuable information gained from
  - Specificity test
  - Stability tests for the prodrug
  - Stability tests for the active drug in the presence of prodrug
- Clinical Samples Handling
  - Sample collection procedure to implement stabilization method if needed
  - Cover the conditions of sample collection and storage

## Consideration during Sample Analysis

- Well controlled and consistent sample analysis procedure
- Analysis within the validated storage conditions
  - Cumulative Benchtop and long term storage
  - F/T cycles
  - Post-extraction storage
- Incurred Sample Reanalysis (ISR)
  - Verify whether the method and stabilization procedure (if applicable) is suitable for a PK study

# Case Example 1



Spiked Concentration of Prodrug /Active drug	% Converted Active Drug	
	Heparin-treated human blood	EDTA-treated human blood
5000/0 ng/mL (ice bath for 1 hour)	25.80%	0.61%
5000/0 ng/mL (ambient temperature for 1 hour)	24.80%	2.88%


## Final Method includes:

- K<sub>3</sub>EDTA as the anticoagulant
- Sample processing at ambient temperature

- Method Validation: QC (1:1 of prodrug/active drug) stability evaluated in plasma, whole blood, and post-extraction solvent
- Sample Analysis:
  - Prodrug/active drug concentration ratio for all study samples < 1:2
  - ISR results
    - Prodrug: 88% passing
    - Active drug: 85% passing

## Treosulfan

## Case Example 2



- Temperature and pH are of crucial importance for stability of the prodrug
- Study Sample drug concentration: Treosulfan: S,S-EBDM at 100-fold
  - 1% ex vivo conversion of treosulfan leads to approximate 100% increase of the SS-EBDM level
- Approach to stabilize treosulfan:
  - Treat human blood samples with 1M citric acid (pH below 5)
  - Minimize clinical sample processing time
  - Sample processing during bioanalysis at 5 C
  - Plasma sample stored at -80 C

<https://doi.org/10.1016/j.jpba.2018.02.049>

Romański, Michał & Głowska, Franciszek. (2018). Journal of Pharmaceutical and Biomedical Analysis. 153. 10.1016/j.jpba.2018.02.049.

# Summary



- Ex vivo conversion of prodrugs to active drugs brings challenges in bioanalysis
  - Assay specificity
  - *Ex vivo* degradation of the prodrugs
- Carefully designed method development and validation tests
  - Addressing key assay challenges early on
  - Avoiding pitfalls caused by ex vivo conversion during sample analysis
  - Details are important to capture issues and evaluate impact

# Challenge Question



1. For bioanalysis, the *ex vivo* conversion from a prodrug to its active drug may occur in sample collection, shipment and storage, but not during sample extraction.

True or False?

**False**

